

Supplemental Table 2. Abs used for flow cytometric phenotyping of CD4+ T cell populations.

CD4+ T Cell Panel					
Marker	Surface/Intracellular Staining^a	Fluorophore	Clone	Vendor	Catalog No.
Live/Dead		Yellow	N/A	Invitrogen	L34959
CD3	S/IC	BV421	SP34-2	BD Biosciences	562877
CD4	S	PerCP-Cy5.5	L200	BD Biosciences	552838
CD8	S	Alexa Fluor 700	RPA-T8	BD Biosciences	557945
CD69	S	PE-Cy7	FN50	BD Biosciences	557745
CCR5	S	PE	3A9	BD Biosciences	550632
Ki-67	IC	FITC	B56	BD Biosciences	556026
TNF- α	IC	BV560	MAB11	BD Biosciences	563418

^aS indicates surface staining, IC indicates intracellular staining.